



MANUFACTURERS OF RAPID DIAGNOSTIC TESTS & REAGENTS

1011 US HWY 22 WEST
PHILLIPSBURG, NJ 08865

Tel.: (908) 213-2900
Fax: (908) 213-2901

**COMPARATIVE EVALUATION OF RAPID HIV TESTS MANUFACTURED WITH
AZOG, INC. TEST COMPONENTS.**

Predicate devices tested were Smart Check, Capillus, Genie II and Determine. Gold standard was EIA and PCR results were also compared.

** This information is reproduced from a telefaxed document produced by: Gede Nigeria, GAIDRI (Gede AIDS and Infectious Disease Research Institute) and dated 24th November 2005.

**PRINCIPAL COLLABORATOR: GEDE AIDS AND INFECTIOUS DISEASES
RESEARCH INSTITUTE GEDE FOUNDATION ABUJA NIGERIA**

SUMMARY

We validated the HIV 1/2 test kits at the GEDE AIDS and Infectious Diseases Research Institute, Abuja Nigeria. The performance compared favorably with the kits that are presently taken as gold standard. It achieved 97.75% sensitivity and a specificity of 100% across a panel of five other test kits. This was above some kits, which are presently on the reference list of our HIV test kit use. On the room temperature performance (Result not included), it achieved 100% compared to Genie II, which could not produce results after 24 hours on the bench. The specifications for individual specimen testing were also achieved. It may be useful in healthcare settings with minimal skill where use of high volume reagents may require additional skill.

We may wish to submit this report for your use.

INTRODUCTION

Enzyme immunoassays (EIA) are the most widely used serologic HIV screening test in developed nations. Fourth generation assays include a wide variety of antigens (HIV-1 group M, HIV-1 group O and HIV-2) and detect IgG as well as IgM, reducing the time between HIV exposure and detection of circulating antibody. The use of EIA in resource-poor settings however is by several factors including the need for well-trained personnel, consistent electric supply, equipment maintenance costs and reagent grade water. This led to the development of rapid, instrument free assay which are now widely used in resource-restrained regions of the world. EIA are still suitable for laboratories performing high-throughput testing. But for smaller labs and "points of service" rapid tests are ideal. There are several HIV rapid test kits now in the

market that are highly sensitive and specific (comparable to that of EIA). Rapid test developed most recently are less reliant on cold chain and provide test results in 15-20 minutes (allowing counseling, testing and return of results in a single visit).

Goal: Validation of HIV rapid test kits to validate its cost effective, sensitivity and specificity and

to explore its possible use in the national algorithm of HIV testing in Nigeria.

Reasons for validation HIV rapid test kits for use in a testing algorithm

Routine formal HIV rapid test kit validation is not routine in Nigeria. Individual HIV rapid tests have been evaluated, but multiple test products have not been evaluated with a single specimen set, which allows for the circulation of sensitivity and specificity of various kits in use. Currently within Nigeria there is no standardized use of HIV rapid tests. A large variety of HIV rapid test products are used throughout the country, not all these are of high quality or suitable for use in Nigeria. More importantly, not all sites conducting HIV testing adhere to a single, universally acceptable procedure for sample collection and testing or proper reagent storage. There are four rapid tests widely used in Nigeria; Smart Check, Capillus, Genie II and Determine. We adapted to evaluate the Rapid HIV test from AZOG, Inc. along with the performance of these kits and reports are now as follows. This is a laboratory based validation providing preliminary results on test performance on a single sample set of well-characterized plasma or serum specimens. A structured laboratory environment with trained, supervised laboratorians is the site of this initial validation. Each individual test was evaluated and sensitivity/specificity calculated. The use of a single sample set of allowed for the circulations of sensitivity and specificity of different proposed tests algorithm.

(Summary is attached)

METHODOLOGY

Standard methodology with laboratory QC and QA, good specimen population distribution and known characteristics, in commercial enzyme immunoassay, (EIA/ELISA, Western Blot (including indeterminate samples) and Polymerase chain reaction (PCR). The population was spread across the four geographical zones (West, North, North Central, North East).

Sample Size: 101 samples were run through the standard steps in kit evaluation procedure, maintaining the specific requirements for each kit.

Selection of HIV tests for validation: This was based on circulating HIV test kits.

Data collection and analysis:

All test results were collected on paper forms and entered into a database (MS Excel) for each kit. The sensitivity and specificity of each rapid test were determined by comparing with positive results with the kit performance. For purposes of validation, the sensitivity of each rapid test was calculated. Specificities were calculated as numbers of test, which are reported as negative by EIA ('true negatives'). Sensitivity and specificity were calculated for each of the test individually and also.

GADRI Evaluation, rapid tests for HIV 1/2

Results of Evaluation of HIV 1/2 Test and other predicate devices
Specimens: Plasma collected in EDTA containing tubes.

Area	Sample ID	W. Blot	PCR	AZOG, Inc	Smart Check	Determine	Capillus	Genie II	EIA
North Central Zone									
1	1609	POS	POS	POS	POS	POS	POS	POS	POS
2	1974 NC	NEG	NEG	NEG	NEG	NEG	POS	NEG	NEG
3	1959	POS	POS	POS	POS	POS	POS	POS	POS
4	2011	POS	POS	POS	POS	POS	POS	POS	POS
5	2091	POS	POS	POS	POS	POS	NEG	POS	POS
6	2231	POS	POS	POS	POS	POS	POS	POS	POS
7	2340	POS	POS	POS	POS	POS	POS	POS	POS
8	2360	POS	POS	POS	POS	POS	POS	POS	POS
9	2399	POS	POS	POS	POS	POS	POS	POS	POS
10	2432	POS	POS	POS	POS	POS	POS	POS	POS
11	2487	POS	POS	POS	POS	POS	POS	POS	POS
12	2465	POS	POS	POS	POS	POS	POS	POS	POS
13	2943	POS	POS	POS	POS	POS	POS	POS	POS
14	2967	POS	POS	POS	POS	POS	POS	POS	POS
15	3067	POS	POS	POS	POS	POS	POS	POS	POS
16	3248	POS	POS	POS	POS	POS	POS	POS	POS
17	3249	POS	POS	POS	POS	POS	POS	POS	POS
18	3908	POS	POS	POS	POS	POS	POS	POS	POS
19	3926	POS	POS	POS	POS	POS	POS	POS	POS
20	4133	POS	POS	POS	POS	POS	POS	POS	POS
Northern Zone									
21	488	POS	POS	POS	POS	POS	POS	POS	POS
22	2050	POS	POS	POS	POS	POS	POS	POS	POS
23	2134	POS	POS	POS	POS	POS	POS	POS	POS
24	2161	POS	POS	POS	POS	POS	POS	POS	POS
25	2368	POS	POS	POS	POS	POS	POS	POS	POS
26	2414	POS	POS	POS	POS	POS	POS	POS	POS
27	2461	POS	POS	POS	POS	POS	POS	POS	POS
28	2642	POS	POS	POS	POS	POS	POS	POS	POS
29	2739	POS	POS	POS	POS	POS	POS	POS	POS
30	2804	POS	POS	POS	POS	POS	POS	POS	POS
31	2836	POS	POS	POS	POS	POS	POS	POS	POS
32	3132	POS	POS	POS	POS	POS	POS	POS	POS
33	3137	POS	POS	POS	POS	POS	POS	POS	POS
34	3248	POS	POS	POS	POS	POS	POS	POS	POS
35	3390	POS	POS	POS	POS	POS	POS	POS	POS
36	3268	POS	POS	POS	POS	POS	POS	POS	POS
37	3566	POS	POS	POS	POS	POS	POS	POS	POS
38	3765	POS	POS	POS	POS	POS	POS	POS	POS
39	3898	POS	POS	POS	POS	POS	POS	POS	POS

Area	Sample ID	W. Blot	PCR	AZOG, Inc.	Smart Check	Determine	Capillus	Genie II	EIA
Western Zone									
40	1972	POS	POS	POS	POS	POS	POS	POS	POS
41	1976	POS	POS	POS	POS	POS	POS	POS	POS
42	1989	POS	POS	POS	POS	POS	POS	POS	POS
43	2020	POS	POS	POS	POS	POS	POS	POS	POS
44	2038	POS	POS	POS	POS	POS	POS	POS	POS
45	2092	POS	POS	POS	POS	POS	POS	POS	POS
46	2203	POS	POS	POS	POS	POS	POS	POS	POS
47	2404	POS	POS	POS	POS	POS	POS	POS	POS
48	2420 NC	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG
49	2511	POS	POS	POS	POS	POS	POS	POS	POS
50	2528	POS	POS	POS	POS	POS	POS	POS	POS
51	2546	POS	POS	POS	POS	POS	POS	POS	POS
52	2630	POS	POS	POS	POS	POS	POS	POS	POS
53	2839	POS	POS	POS	POS	POS	POS	POS	POS
54	2938	POS	POS	POS	POS	POS	POS	POS	POS
55	3139	POS	POS	POS	POS	POS	POS	POS	POS
56	3165	POS	POS	POS	POS	POS	POS	POS	POS
57	3321	POS	POS	POS	POS	POS	POS	POS	POS
58	3777	POS	POS	POS	POS	POS	POS	POS	POS
59	3912	POS	POS	POS	POS	POS	POS	POS	POS
Northeastern Zone									
60	2	POS	POS	POS	POS	POS	POS	POS	POS
61	3	POS	POS	POS	POS	POS	POS	POS	POS
62	5	POS	POS	POS	POS	POS	POS	POS	POS
63	6	POS	POS	POS	POS	POS	POS	POS	POS
64	15	POS	POS	POS	POS	POS	POS	POS	POS
65	18	POS	POS	POS	POS	POS	POS	POS	POS
66	20	POS	POS	POS	POS	POS	POS	POS	POS
67	21	POS	POS	POS	POS	POS	POS	POS	POS
68	22	POS	POS	POS	POS	POS	POS	POS	POS
69	27	POS	POS	POS	POS	POS	POS	POS	POS
70	168	POS	POS	POS	POS	POS	POS	POS	POS
71	1475	POS	POS	POS	POS	POS	POS	POS	POS
72	2053	POS	POS	POS	POS	POS	POS	POS	POS
73	2385	POS	POS	POS	POS	POS	POS	POS	POS
74	2549	POS	POS	POS	POS	POS	POS	POS	POS
75	2855	NEG	POS	POS	POS	POS	POS	POS	POS
76	3126	NEG	POS	POS	POS	POS	POS	POS	POS
77	3213	NEG	POS	POS	POS	POS	POS	POS	POS
78	2865	NEG	POS	POS	POS	POS	POS	POS	POS
79	3909	NEG	POS	POS	POS	POS	POS	POS	POS
80	4207	NEG	POS	POS	POS	POS	POS	POS	POS
81	2420 NC	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG
82	3685 NC	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG
83	23 NC	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG
84	13 NC	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG

Area	Sample ID	W. Blot	PCR	AZOG, Inc.	Smart Check	Determine	Capillus	Genie II	EIA
85	3314 NC	POS	NEG	NEG	NEG	POS	NEG	NEG	NEG
86	3660 NC	POS	NEG	NEG	NEG	NEG	NEG	NEG	NEG
87	3760 NC	POS	NEG	NEG	NEG	NEG	NEG	NEG	NEG
88	3768 NC	POS	NEG	NEG	NEG	NEG	NEG	NEG	NEG
89	3837 NC	POS	NEG	NEG	NEG	NEG	NEG	NEG	NEG
90	4239 NC	POS	NEG	NEG	NEG	NEG	NEG	NEG	NEG
91	3394 ID	Indeterminate	POS	POS	NEG	NEG	NEG	NEG	NEG
92	3975 ID	Indeterminate	POS	NEG	NEG	NEG	NEG	NEG	NEG
93	2981 ID	Indeterminate	POS	NEG	POS	POS	POS	POS	POS
94	3210 ID	Indeterminate	POS	POS	POS	POS	POS	POS	POS
95	3230 ID	Indeterminate	POS	POS	POS	POS	POS	POS	POS
96	3763 ID	Indeterminate	POS	POS	POS	POS	POS	POS	POS
97	2428 ID	Indeterminate	POS	POS	POS	POS	POS	POS	POS
98	2485 ID	Indeterminate	POS	POS	POS	POS	POS	POS	POS
99	2468 ID	Indeterminate	POS	POS	POS	POS	POS	POS	POS
100	3243 ID	Indeterminate	POS	POS	POS	POS	POS	POS	POS
101	4230 ID	Indeterminate	POS	POS	POS	POS	POS	POS	POS

Results of the evaluation of HIV test kits using test components available from AZOG, Inc. and other predicate devices.

Panel System	Sensitivity TP/TF + FNEG x 100	Specificity TNEG/TNEG + FPOS x 100	Positive Predictive Value N. TPN/N. TP + NFPOS x 100
AZOG, Inc HIV – 1/2 Kit	97.75%	100.00%	100.00%
Smart Check HIV – 1/2 Kit	97.75%	100.00%	100.00%
Determine HIV – 1/2 Kit	97.75%	92.31%	96.96%
Capillus HIV – 1/2 Kit	97.65%	100.00%	100.00%
Genie II HIV – 1/2 Kit	97.75%	100.00%	100.00%
EIA HIV – 1/2 Kit	97.75%	100.00%	100.00%

REFERENCE STANDARD

**TRUE POSITIVE = 89, TRUE NEGATIVE = 12
TOTAL = 101**

AZOG, Inc. HIV – 1/2 Kit	True positive = 87, True negative = 12; False positive = 0, False negative = 2
Smart Check HIV – 1/2 Kit	True positive = 87, True negative = 12; False positive = 0, False negative = 2
Determine HIV – 1/2 Kit	True positive = 87, True negative = 12; False positive = 1, False negative = 2
Capillus HIV – 1/2 Kit	True positive = 87, True negative = 12; False positive = 0, False negative = 3
Genie II HIV – 1/2 Kit	True positive = 87, True negative = 12; False positive = 0, False negative = 2
EIA HIV – 1/2 Kit	True positive = 87, True negative = 12; False positive = 0, False negative = 2